

QUAT KIT



Code 7057-01 | Drop Count, 1 Drop = 2, 5, 10 ppm

QUANTITY	CONTENTS	CODE
15 mL	*Phenolphthalein Indicator, 0.5%	*2258-E
15 mL	*Sulfuric Acid, 0.5 N	*6090-E
60 mL	Quat Titrating Solution	3996-H
30 mL	Toluidine Blue O Indicator	3995-G
60 mL	*EDTA Solution	*7117-H
1	Test Tube, 5-10-15-20-25-30 mL, plastic, w/cap	0715
1	Quat/Polyquat Endpoint Color Chart	3613-CC

To order individual reagents or test kit components, use the specified code number.




*Reagent is a potential health hazard. **READ SDS:** lamotte.com
Emergency information:
Chem-Tel USA 1-800-255-3924
Int'l, call collect, 813-248-0585



PROCEDURE

1. Rinse test tube [0715] with sample water. Fill with desired sample size selected from the table.
- | Sample Size | Equivalence (ppm Per Drop) |
|-------------|----------------------------|
| 25 mL | 1 drop = 2 ppm |
| 10 mL | 1 drop = 5 ppm |
| 5 mL | 1 drop = 10 ppm |
2. Add 5 drops *EDTA Solution [7117]. Swirl to mix. NOTE: If the hardness of the sample is greater than 500 ppm, add 5 more drops of EDTA Solution.
3. Add 2 drops of *Phenolphthalein Indicator, 0.5% [2258]. Swirl to mix. If colorless, proceed to Step 4. If pink, add *Sulfuric Acid, 0.5N [6090] dropwise, until the pink color disappears.
4. Add Toluidine Blue O Indicator [3995] as follows:
- | | |
|--------------|-------------|
| 25 mL sample | add 8 drops |
| 10 mL sample | add 3 drops |
| 5 mL sample | add 2 drops |
- Swirl to mix. Sample should turn light blue.
5. While swirling test tube, add Quat Titrating Solution [3996] one drop at a time, until color changes from blue to purple. Hold bottle vertically. For best

results, when the color change is first detected, use the Endpoint Color Chart [3613-CC] as shown to match the color of the solution exactly to the endpoint. Continue adding Quat Titrating solution one drop at a time until color matches endpoint. Count the number of drops added.



6. Calculate result in ppm.
- 25 mL sample:** subtract three [3] from number of drops used in Step 5. Multiply by 2.
- 10 mL sample:** subtract one [1] from number drops used in Step 5. Multiply by 5.
- 5 mL sample:** subtract one [1] from number drops used in step 5. Multiply by 10. Record as ppm Quat.
- NOTE: The quat equivalence is based on n-alkyldimethylbenzyl ammonium chloride, molecular weight 360. If a quat of different molecular weight is tested, multiply the equivalence by:

$$\text{equivalence} \times \frac{\text{molecular weight}}{360}$$