



5-WAY Nitrate, Nitrite, pH, Alkalinity and Total Hardness

25 Test Strips

Code 3038-G

Insta-Test[®] ANALYTIC

Immerse strip. SWIRL 2 TIMES.
Remove with pads face up.
DO NOT SHAKE OFF EXCESS WATER. Immediately read **Hardness > Alkalinity.**
Wait 30 seconds, read **pH > Nitrite > Nitrate.**

- Keep wet fingers out of the vial
- Close vial tightly after removing strip
- Store in a cool, dry place

LaMotte Company
Chesterstown, Maryland 21620 USA
www.lamotte.com

NO₃
Nitrate
ppm

Fresh & Salt	0	20	40	OK	80	160	200

NO₂
Nitrite
ppm

Fresh & Salt	0	0.5	1	3	5	10

pH

Fresh	6.0	6.5	7.0	OK	7.5	8.0	8.5	9.0
Salt				OK				
	6.0	6.5	7.0	7.5	8.0	8.5	9.0	

Alk
Alkalinity
ppm

Fresh & Salt	0	40	80	OK	120	180	240

TH
Hardness
ppm

Fresh	0	30	60	OK	120	180



8

46306

37335

6

Drop holding box and dotted tolerance line